SMART  TOOLS

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Course Description

Smart Tools Plus (STP) education and tools are a form of instrument-assisted soft tissue mobilization (IASTM) therapy. The five stainless steel instruments were designed for the unique programming taught in our classes. The purpose of this Level 1 training is to introduce clinicians to STP programming and tools and workshop each and every region of the body. Level 1 training includes assessment, treatment, physiological effects/benefits and indication/contraindications to IASTM.

Level 1 Course Objectives

The primary objective of Level 1 is to develop an understanding of STP and how to apply it into the full spectrum of musculoskeletal treatment approaches. By the end of the training, students will be able to:

- Demonstrate a working knowledge of the STP instruments, STP treatments and potential effects and benefits.
- Identify and discuss the indications, contraindications (relative and absolute) of IASTM.
- Review and develop a better understanding of soft tissue injury, healing and potential reactivity to instrument-assisted soft tissue mobilization (therapeutic and adverse).
- Develop skill and competence in the STP application of IASTM to the major regions of the spine and extremities.

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Level 1 Training Content

Hour 1: Introduction to the Technique

History, physiological rationale of IASTM, principles relating to soft tissue dysfunction and anticipated results from using STP instruments & Programming are covered. A review of the current scientific literature and evidence based model is discussed.

Hours 2 – 3: Introduction of STP Instruments and Treatment Philosophy

Demonstration of each instrument and edge, proper handholds and treatment strokes are taught. Lab instruction includes 'scanning', 'pinning', 'prying & lifting', and 'clearing the path'.

Hours 4 – 8: Lower Extremity: Ankle/Foot, Lower Leg (anterior, posterior, deep posterior, and lateral), Knee

Participants learn how to apply the STP treatment strokes to the Lower Extremity: Below the knee. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

Upper Leg: Anterior, Posterior, Lateral, and medial thigh.

Participants learn how to apply the STP treatment strokes to the Lower Extremity: Upper Leg. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

Hip:

Participants learn how to apply the STP treatment strokes to the hip external rotators and tensor fascia latae. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

Hours 8 – 16: Torso: Lower Back (Iliac Crest, Sacroiliac joint, Paraspinals, Quadratus Lumborum

Participants learn how to apply the STP basic treatment strokes to the lower back. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region.

Participants practice the above in lab format.

Torso: Upper Back, Chest, Abdomen

Participants learn how to apply the STP treatment strokes to the rhomboids and latissimus dorsi. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

Torso: Chest/Abdomen

Participants learn how to apply the STP treatment strokes to the Pectoralis and Abdomen region. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

Head and Neck.

Participants learn how to apply the STP treatment strokes to the head and neck region including cervical paraspinals, upper trapezius, levator scapulae, and sternocleidomastoid. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

Upper Extremity: Shoulder, Elbow and Wrist:

Participants learn how to apply the STP treatment strokes to the upper extremity. Suggested patient and clinician positioning and instrument choice is demonstrated and practiced for this region. Participants practice the above in lab format.

TREATMENT EFFECTS OF IASTM

Break up abnormal densities in tissue

Stimulate nervous system

Reinitiate first-stage healing in the body via proliferation

Tool Assisted Soft Tissue Technique White Paper (1/2)

INTRODUCTION

Optimal human movement requires a careful balance of both mobility and stability of the joints and regions to disseminate forces throughout the body and avoid injury. Grey Cook, PT describes a joint-by-joint philosophy of how some areas of the body primarily require mobility and some areas of the body primarily require stability. {Cook, 2010 #377} Figures 1 & 2 are a graphical representation of the joint by joint approach.

When a region or joint that requires mobility has less than adequate joint range of motion (ROM), the restriction may be caused by either soft tissue tightness or joint restriction. Joint restriction can be addressed with manual therapy techniques like joint mobilization and manipulation. In addition, home exercises can be prescribed with banded joint mobility to improve joint mechanics. {Starrett, 2013 #378}

Soft tissue restrictions can be addressed in many different ways. Some of those ways are massage, dry needling, pin and stretch, therapeutic ultrasound, heat, and stretching to name a few. A unique way for providers to address soft tissue restrictions is Instrument Assisted Soft Tissue Manipulation (IASTM).

The purpose of this white paper is to discuss IASTM performed by clinicians in order to improve range of motion by removing soft tissue adhesion and scar formation, improve motor control and stimulate mechanoreceptors that lay underneath the surface of the skin. Clinicians include massage therapists, physical therapists, athletic trainers, chiropractors, and physicians among other providers who are licensed to provide hands on manual therapy techniques in which this technique falls within their scope of practice.

THE CHALLENGE:

A patient presents to the clinician's office complaining of pain that is causing some restriction in activity or work or causing disruption to the patient's life. The clinician attempts to provide a diagnosis based on the history, mechanism of injury, and medical history. The evaluation of the injury includes range of motion testing. A good clinician will attempt to figure out what is causing the restriction in range of motion. One common cause of restricted range of motion is inadequate tissue flexibility.

- Inadequate tissue flexibility increases susceptibility to both overuse syndromes and acute injuries. {Sainz de Baranda, 2010 #366}
- Improvement of flexibility is often a goal of interventions for injury prevention, performance enhancement, and injury rehabilitation. {Sainz de Baranda, 2010 #366}
- A variety of stretching techniques and heating modalities are used routinely to promote flexibility. A theorized cause of apparent tissue tightness is tissue extensibility dysfunction (TED). {Cook, 2010 #377}
- Musculoskeletal injury, even if microtraumatic, can produce scar tissue or adhesions that alters the properties of collagenous tissues. The term neurodynamics refers to peripheral nerve sliding and/or tension development that normally occurs without symptoms. Appropriate neurodynamic treatment strategies (i.e., "slides" or "tensioners") may alleviate neural symptoms that arise from TED. {Butler, 2000 #379}
- Neuromuscular activation patterns may involve strong contractions that produce joint compression and excessive muscle tension. {Cook, 2010 #377}
- When treating TED, clinicians need to identify its cause in order to implement an intervention that will be effective in producing long-term flexibility improvement. IASTM is a therapeutic technique that is based on the soft tissue mobilization rationale introduced by James Cyriax. IASTM differs from traditional cross-friction or transverse friction massage. Specially designed instruments are used to apply longitudinal pressure along the course of muscle fibers, and treatment typically includes application to more than the tissues at the isolated location of pain. The instruments are thought to facilitate the clinician's ability to detect altered tissue properties, as well as facilitate the patient's awareness of altered sensations within the treated tissues. Increased vibration within the instrument is believed to be an indication of abnormal tissue properties. Additionally, the instruments are believed to provide a mechanical advantage that allows the clinician to achieve greater depth of mechanical force transmission than that which can be produced with the hands, while also reducing compressive stress on the clinician's hands. The clinical use of soft tissue mobilization instruments is purported to enhance treatment effectiveness, {Aspegren, 2007 #360; Bayliss, 2011 #356; Burke, 2007 #361; Davidson, 1997 #375; Fowler, 2000 #376; Hammer, 2005 #362; Howitt, 2009 #358; Laudner, 2014 #354; Lee, 2014 #353; Loghmani, 2009 #349} particularly in areas of fibrosis.
- The inducement of tissue microtrauma is believed to elicit a local inflammatory response that promotes breakdown of scar tissue, release of adhesions, synthesis of new collagen, and connective tissue remodeling. {Hreljac, 2000 #367}
- IASTM treatment of enzyme-induced tendinitis in rats has been shown to promote fibroblast proliferation, collagen synthesis, collagen maturation, and collagen alignment. {Davidson, 1997 #375}
- Functional benefits included increased stride length and stride frequency. {Bayliss, 2011 #356}
- A study of ligament healing in rats following IASTM treatment documented greater cellularity, improved collagen alignment, fewer adhesions and granular tissue, and increased ligament strength and stiffness. {Loghmani, 2009 #349}

MOBILITY PATTERNS


Figure 1: Joint/Regions That Primarily Require Mobility
Big Toe, Ankle, Hip, Thoracic Spine, Glenohumeral Joint, Upper Cervical Spine (C0, C1, C2)

MOTOR CONTROL PATTERNS


Figure 2: Joints/Regions That Primarily Require Stability
Arch, Knee, Low Back, Scapulae, Lower Cervical Spine, Elbow/Wrist

- *Clinical studies of IASTM administered to patients with tendinopathies have demonstrated pain resolution, improved ROM, and return to normal function at a faster rate than that observed for natural healing and traditional therapeutic interventions. Case studies have documented similar IASTM benefits, but other therapeutic interventions also were administered in those cases.*{Bayliss, 2011 #356;Burke, 2007 #361;Hammer, 2005 #362;Howitt, 2009 #358}

The goal of therapy is to provide an optimal environment for healing, by either modifying physiologic responses to injury (e.g., inflammation, muscle spasms, pain) or enhancing components of normal musculoskeletal function (e.g., increase range of motion, increased muscular strength).{Starkey, 2004 #380} *No single therapeutic approach in isolation will completely resolve an impairment of musculoskeletal function. Smart Tool Plus® represents an approach to soft tissue manipulation that uses five different stainless steel instruments to release scar tissue, adhesions, and fascial restrictions. Smart Tools Plus is designed to be used in conjunction with a thorough medical evaluation, movement based examination and other modalities to help patients overcome injuries.*

Therapeutic stretches and exercises (preferably eccentric exercises){Nelson, 2004 #372} are necessary to promote tissue lengthening and collagen fiber realignment, which help to prevent the released tissue from becoming restricted again.

Table 1 outlines indications for use of ISTM, which is primarily focused on soft tissues; Currently, the benefits derived from ISTM include release of fascial restrictions, breakdown of collagen cross-linkages, increased blood flow, and possibly an increase in regenerative cellular activity.{Davidson, 1997 #375;Nelson, 2004 #372;Portillo-Soto, 2014 #363;Vardiman, 2015 #352} *Even with the greatest care, however, there are potentially adverse treatment responses to ISTM that may occur. The patient may experience discomfort during administration of the treatment, and petechiae (i.e., bruising) may become apparent during or after the treatment (Figure 2). Bruising results from localized trauma, which may be associated with separation of adhesions from healthy tissue.*{Hammer, 2008 #359}

The patient must be informed about the potential effects and benefits of the treatment, and other therapeutic modalities should be administered for pain management. See Table 2 for precautions and contraindications.

The SMART Tools® have the desired weight, sharpness, and ergonomic design that practitioners prefer for optimal results. The edges are sharper than most IASTM tools in order to maintain tissue tension throughout the treatment strokes and to decrease overall treatment time by creating the desired localized inflammatory process needed for tissue remodeling/healing to take place. There are also single and double beveled edges on various tools. The single beveled edge allows for more tissue tension and penetration, while the double beveled edge allows for more efficient treatment sessions due to allowing the practitioner to perform treatment strokes in a variety of directions.

Therapeutic Outcome:

Research evidence that documents the effectiveness of ISTM is limited; however, several experimental and case studies have produced positive findings. McLaughlin{McLaughlin, 2006 #462} investigated the effectiveness of ISTM for reduction of edema associated with ankle sprains. Participants were intercollegiate, intramural, and high school athletes. They were randomly assigned to either a traditional edema control protocol or a traditional edema control protocol that was combined with ISTM. There was no significant difference in edema control between the two groups (i.e., both therapeutic protocols were equally effective in reducing edema). Achievement of full weight-bearing status for the group that received ISTM averaged one day sooner than that for the comparison group, however. Lohmani and Warden{MT, 2009 #463} assessed short-term and long-term effects of ISTM on healing of the medial collateral ligament (MCL) in an animal model. Fifty-one rats had bilateral MCL injuries surgically induced to both hind legs. Thirty-one of the animals received instrument-assisted cross-fiber massage treatment three times per week for three weeks to one extremity (i.e., the other was untreated), and 20 animals received the same treatment for 10 weeks to one extremity. At four and 12 weeks post injury, MCLs were harvested for testing. The non-treated MCLs were noted to have more adhesions and granular tissue, which made the harvesting process more challenging. In comparison to non-treated MCLs at four weeks post injury, the treated MCLs had 43.1% greater tensile strength, 39.7% greater stiffness, and were able to absorb 57.1% more energy before failure. At 12 weeks post-injury, the treated MCLs were 15.4% stiffer than the non-treated ligaments, but there was no significant difference in tensile strength or energy absorption to failure. The treated ligaments demonstrated greater cellularity and better collagen fiber alignment when compared to the nontreated ligaments at weeks 4 and 12. The non-treated ligaments demonstrated greater scarring and more poorly organized collagen, especially at four weeks postinjury.

Use of ISTM was reported in a case study that involved a 59-year-old man with a one-year history of intense low back pain (i.e., subacute lumbar compartment syndrome) that caused him to miss two to three days of work every two to three months.{Hammer, 2005 #362} *His initial treatment protocol consisted of bed rest and analgesics. His pain was managed well enough to allow him to work as a shoe salesman until he experienced pain that would not subside, and that prevented performance of his activities of daily living for a period of two weeks. ISTM treatments were administered to the hamstrings, sacrum, right hip lateral rotators, and low back region. The patient received six treatment sessions (twice per week for three weeks), and each included performance of two sets of three stretches to the affected area after administration of ISTM. The patient was instructed to perform the stretches at home between treatment sessions. After the six sessions, the patient was asymptomatic and able to complete all tasks necessary for daily living and work.*

CONCLUSION

The SMART Tools were developed to assist in effective treatment of adhesions and scarring brought on by surgeries, injuries, and overuse of muscles. Our tools will assist Medical Doctors, Chiropractors, Acupuncturists, Physical Therapists, Occupational Therapists, Massage Therapists, Athletic Trainers and others, to effectively treat patients with less fatigue to the practitioner. We have experimented with various tools and confident the SMART Tissue Tools will provide years of satisfaction at a fraction of the cost of comparable products on the market. Our engineered double edge technology applied through the use of precision machining creates an unrivaled contact with the skin surface and the tool to transmit a unique vibration which helps isolate the troubled area. The medical grade 304 Stainless Steel, polished to a mirror-like finish, allows the use of topical solutions to be used without any risk or damage the tools. This will enable the practitioner to create a controlled inflammatory process which will allow the tissue to begin the healing process, breaking down scar tissue and adhesions in the troubled area and accelerate the recovery process.

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